


1. Introduction

The Distortionless Enhancement by Polarization Transfer (DEPT) sequence uses polarization transfer from protons to other nuclei via one covalent bond to increase signal strength. The experiment is typically used for CH_n multiplicity determination. The ¹³C DEPT-135 experiment on cholesterol is described herein. This example demonstrates the basic procedure of double resonance 1D NMR data acquisition and processing on Tecmag spectrometers.

2. Pulse sequence



WALTZ-16 = RR $\bar{R}\bar{R}$ \bar{R} RR \bar{R} RR $\bar{R}\bar{R}$ \bar{R} RRR,
 R = 90_x180_{-x}270_x

(b)

Pulse width and phase cycle:

P1 (H90°): 0

P2 (H180°): phH2 = 0, 2, 1, 3.

P3 (H135°): phH3 = 1, 1, 1, 1, 3, 3, 3, 3.

P4 (C90°): phC1 = (0)₈, (2)₈, (1)₈, (3)₈.

P5 (C180°): phC2 = (0, 2)₄, (1, 3)₄.

Receiver: phRX = (1)₂, (3)₄, (1)₂, (2)₂, (1)₄, (2)₂, (3)₂,
 (1)₄, (3)₂, (0)₂, (2)₄, (0)₂.

(phH2, phH3, phC1, phC2, phRX are 1D phase tables.

All tables are in 4 step mode.)

Event Number	1	2	3	4	5
Name:	R1	R2	R3	R1	R2
Delay	90	180	270	90	180
F1_Ampl					
F1_PhMod	X	-X	X	X	-X





Fig. 1. (a) The ¹³C DEPT-135 sequence with WALTZ-16 sequence for ¹H decoupling. (b) The sequence realized in the NTNMR sequence editor.

3. Experiment


Sample:	Cholesterol in CDCl ₃ (50mg/ml)
Spectrometer:	7 Tesla Magnet with Tecmag HF3 discovery
Probe:	Nalocac D300-5 OWB 5mm ¹ H/ ¹³ C Switchable probe
¹ H hard pulse:	13.9 kHz (90° = 18 μs @ 5W)
¹ H decoupling field:	5.6 kHz (90° = 45 μs @ 800 mW)
¹³ C hard pulse:	59.5 kHz (90° = 4.2 μs @ 250 W)
τ:	3.5 ms (= 1/2J _{C,H} = 140 Hz)
SW +/-:	± 6.5kHz
Last Delay:	2s
Scans 1D:	512



Notes:

1. Before editing the sequence (Fig. 1b), calibrate the 90° pulse widths of ¹H and ¹³C using the nutation experiment (see note, "One Pulse Experiment and Pulse Calibration").
2. Set up the WALTZ sequence according to the note, "¹³C NMR Spectra with ¹H WALTZ Decoupling".
3. The center of pulses P2 and P4 (also P3 and P5) should be aligned. Since P2 > P4 (and P3 > P5) P2 (and P3) have to split into 3 pulses. The delay of P2's (and P3's) middle pulse equals to P4 (and P5), and the delay of both sides is (P2 - P4)/2 [and (P3 - P5)/2]. The middle pulse of P2 (and P3) falls on the same event as P4 (and P5).

4. Results



CH Fig. 2. The ¹³C DEPT-135 spectrum of cholesterol obtained using the sequence shown in Fig. 1. Since signals are generated from proton polarization, both quaternary carbons and the solvent peak do not appear.

5. References

1. M.R. Bendall, D.M. Doddrell, D.T.Pegg, *J. Am. Chem. Soc.* **1981**, 103, 4603-4605.
2. D.M. Doddrell, D.T.Pegg, M.R. Bendall, *J. Magn. Res.* **1982**, 48, 323-327.
3. S. Braun, H.-O. Kalinowski, S. Berger, "150 and More Basic NMR Experiments", Wiley-VCH, 1999, p.180-182.