

Protasis High Throughput NMR Structure Confirmation System

Typical Configuration

- ✓ Protasis CapNMR™ ICG Probe
- ✓ Protasis One-Minute NMR™ Software
 - Talks to Varian or Bruker Spectrometers
 - Seamless Interface to ACD Structure Confirmation Software
 - Accurate %Yield and Purity in a single step with ERETIC
- ✓ Protasis MultiSolvent HTSL + DSCM
- ✓ CTC/LEAP PAL Liquid Handler
- ✓ Vacuum Centrifuge
- ✓ Varian or Bruker Spectrometer
- ✓ 400 MHz Shielded Magnet

